



MDS SERIES DISTILLING DEVICE

*Production of
Custom Sizes and
Features

MDS Series Water Distillers are devices that help to turn mains water into pure water in the fastest and ergonomic way by passing it through certain processes and are used in many are

CONVENIENT USE

- It has a stainless steel heater,
- All surfaces in contact with water and steam are made of stainless material.
- It can be used as table top or wall mounted type.
- It consists of an ultra robust design designed to last for many years.

RELIABLE DESIGN

- Empty boiling kettle and protection against water shortage, Low water level warning,
- High efficiency even at low water pressure,
- Easy and understandable operation,
- Illuminated on/off switch indicating that the device is working
- Gas safety thermostat against overheating.



	MDS-4	MDS-8	MDS-12
CAPACITY	4 Lt / Saat	8 Lt / Saat	12 Lt / Saat
HEATER	3000 W	6000 W	9000 W
COOLING WATER	40 Lt	80 Lt	120 Lt
OUTER DIMENSIONS (EXBXD) mm	400x622x400	400x622x400	400x622x400
MATERIAL	INTERIOR MATERIAL	Stainless Steel	
	EXTERIOR MATERIAL	Electrostatic Powder Coated Steel Against Rust	
SAFETY	Over Temperature Protection		
	Water Level Sensor		
NET WEIGHT (kg.)	26	26	26
PACKED WEIGHT (kg.)	30	30	30

OPTIONAL REQUIREMENTS

- Pre-treatment possible,
- Effluent conductivity 2-5 µm.

