



MDK SERIES DESICCATOR CABINETS

*Production of
Custom Sizes and
Features

VARIOUS SIZES

MDK Series Desiccator cabinets can be produced with the desired width, length, depth and number of shelves according to user needs.

ROBUST DURABILITY

- Protects the shelf life, quality and durability of sensitive products.
- Provides a protected area for precise measurements and testing.
- It ensures the protection of materials from the benefits arising from moisture.
- Reduces costs by preventing product loss.

OPTIONAL REQUESTS

- Production in desired dimensions and sizes,
- Extra shelf and sample container.



	MDK-50
EXTERIOR DIMENSIONS (WxLxD)	325x420x580 mm
INTERIOR DIMENSIONS (WxLxD)	280x50x520 mm
VOLUME	50 liter
APPROXIMATE WEIGHT	8 kg
MAXIMUM WEIGHT THAT THE SHELVES CAN CARRY	6 kg
CABIN MATERIAL	Acrylic Sheet
FRAME MATERIAL	Aluminium, Polymethylmethacrylate Based Thermoplastic Material
STANDARD SHELF COUNT	4
INDICATOR	Digital Temperature and Humidity
OPERATING TEMPERATURE	-20 / + 70 °C
USAGE MANUAL	Room Conditions
OTHER FEATURES	Bottom Collection Tray, Silica Gel

