



# MIPROLAB

## MKP SERIES NITROGEN PROTEIN DISTILLATION UNIT

\*Production of  
Custom Sizes and  
Features

### ENHANCED CONTROL SYSTEM

- ▶ The device body is completely resistant to all sweat agents and does not cause any spillage.
- ▶ The device is programmable microprocess controlled.
- ▶ The display of the device is LCD, the display is alphanumeric and the keys are touch buttons.
- ▶ At the end of the distillate, waste initiation is activated automatically.
- ▶ The device warns the user by giving warning messages.
- ▶ The device notifies the user if there is water in the steam generator.
- ▶ The amount of alkali to be added can be adjusted.
- ▶ The steam generator has a power of 1700 Watt.
- ▶ The device is complete with accessories.



### MKP

DIMENSIONS (EXBXD) mm

505x713x355

NET WEIGHT (kg.)

36

PACKED WEIGHT (kg.)

40

**MIPROLAB**  
MIPROLAB LABORATORY EQUIPMENT INC.

AĞAÇ İŞLERİ SANAYİ SİTESİ 1362. CAD. NO:41 ANKARA/TÜRKİYE

Phone: +90 312 395 75 55 | Whatsapp: +90 542 395 75 56

E-mail: info@miprolab.com.tr

