



MT SERIES TIK TIK DEVICE

*Production of
Custom Sizes and
Features

MT Series TIK TIK devices are designed to compact all kinds of soil samples easily and quickly.

POWERFUL SYSTEM

Table top and single person portable crate type.
Suitable for compaction of soil sample.
5 sample containers are supplied with the device.
Convenient use with manual numberer.
220 V, 50 Hz



MT	
DIMENSIONS (EXBXD) mm	300x345x500
NET WEIGHT (kg.)	20
PACKED WEIGHT (kg.)	23

STABILITY AGGREGATE DEVICE

Convenient use with digital display.
Table top type and can be carried by one person.
Suitable for stabilising 6-8-12 soil samples.
Adjustable number of dives and automatic stop when the number of dives is completed.
Possibility of use with 0.060 - 0.125 - 0.25 - 0.5 - 1 - 2 - 4 and 5 mm sieves. Ability to attach and detach one sieve series to another sieve series.
Water reservoir large enough to submerge 4 sieve series and 6-8-12 water reservoirs to submerge the sieve set.

