



MTR SERIES THERMORACTOR

*Production of
Custom Sizes and
Features

OPTIONAL REQUESTS

- Possibility to operate at different temperatures between 5 °C and 200 °C,
- 30, 60, 90, 120 minutes time setting or continuous operation,
- Automatic switch-off and audible warning at the end of the reaction time,
- Possibility to work with 13-16-22 mm diameter glass tubes, 6-8-12-24 tubs,
- Optional separately controllable double block,
- Safety control against overheating,
- Changeable block,
- Block suitable for different tubes,
- Touch screen option.



	MTR TEKLİ	MTR ÇİFTLİ
INSTALLED POWER	500 W	1000 W
DIMENSIONS (EXBXD) mm	220x135x305	220x135x440
NET WEIGHT (kg.)	5	10
PACKED WEIGHT (kg.)	7	12

