



MIPROLAB

MUC SERIES ULTRASONIC WATER BATH

**Production of
Custom Sizes and
Features*

MUC Series ultrasonic baths are used for cleaning (technical glassware such as burettes, pipettes, petri dishes), medical instruments, all kinds of metal parts, electronic parts, degassing of liquids for concentration determination, acceleration of dilution processes, simultaneous cleaning and disinfection, emulsion preparation, sample preparation for analysis and many other areas.

WIDE FEATURE AREA

Provides deep cleaning of glass materials, metal products, electronic components and various apparatus, Robust one-piece vibration tank and body made of stainless steel, Continuous and effective ultrasonic wavelength guarantees precise and precise cleaning on all surfaces, Control unit provides temperature and time control. Displays set and actual values in parallel, In case of frequently repeated operation, the memory function stores the stored values and activates them immediately on demand.


STRONG SYSTEM

- The heating system ensures that the desired temperature is reached in a very short time,
- The temperature inside the bath can be adjusted between ambient temperature and 80 °C,
- Time control can be set between 1-30 minutes.

MUC-45	
DIMENSIONS (EXBXD) mm	530x470x325
NET WEIGHT (kg.)	20
PACKED WEIGHT (kg.)	23

OPTIONAL REQUESTS

- Production in Special Dimensions and Specifications (Dimension/Software/Programme),
- 4.3 / 5.5 and 7.2 Inch Screen Options.
- Internal Memory (125 Days),
- Optional accessories.
- Internal and external circulation addition.



MIPROLAB

MIPROLAB LABORATORY EQUIPMENT INC.

AĞAÇ İŞLERİ SANAYİ SİTESİ 1362. CAD. NO:41 ANKARA/TÜRKİYE

Phone: +90 312 395 75 55 | Whatsapp: +90 542 395 75 56

E-mail: info@miprolab.com.tr

